QuickLase®

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QuickLase® 810 nm

All the QuickLase systems are small, lightweight and portable.
Battery (up to 4 hours of use) and mains powered.
Ergonomic system for comfortable easy use.
Touch screen with presets, multi-language menus and very easy to use.

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Cutting edge technology at an easily affordable price.

QuickLase offers the ability to cut, coagulate and sterilise. It reduces both tissue scarring and post-operative discomfort. Healing is rapid and more aesthetic. It provides the perfect outcome for patients and dentist alike.

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Curretage / Curetaje

Soft Tissue system

Points to consider before buying a laser:

1. Wavelength

The wavelength used is the most important characteristic. It determines how the laser light interacts with the soft tissue. Notably, the absorption by the appropriate tissue chromophores and the penetration depth into the tissue. When using only one wavelength of laser light, 810 nm is one of the most versatile overall. QuickLase always uses 810 nm in all our models as it has a greater effect on coagulation whilst cutting soft tissue.

2. Maximum Power

The output power determines the speed with which the various procedures can be completed.

In dentistry, the hand speed and skill of the practitioner will determine the required power output for that procedure. A higher power will allow for the treatment to be completed faster than a lower power but all procedures can be successfully completed with all our lasers. QuickLase offers you both 3 & 5 watt (810nm) lasers, and 10 & 20 watt (810 + 980nm) systems.

3. Emission modes / Usage

The laser can be set to either Continuous Wave or Pulsed Mode (modulated), this relates to the type of output. The prime benefit of the pulsed mode is the capacity of the tissue to cool between successive pulses. This allows for a reduction in the amount of, or in some cases the need for anaesthesia required. This mode is however slower than Continuous Wave mode. Before choosing an emission mode you must also consider the power setting, as the average power with the continuous mode is the power selected, but when you set the pulse mode the average power will be less. The difference is determined by the pulse frequency. That's why you must allow more time when using the pulse mode than continuous mode at the same power output. QuickLase has continuous and pulse modes.

4. Delivery medium

Most dental diode lasers have a flexible fibre optic to deliver the laser output. There are a number of things to consider when using a fibre optic, but the most important is the diameter of the fibre. Using a smaller diameter fibre will increase the power density at the fibre tip, so a smaller diameter fibre will require less power than a larger diameter, and a larger diameter fibre will cover a bigger area. QuickLase offers two fibre optic diameters, 200 µm and 400µm (microns).

5. Price

Of course a determining factor in making a decision is the price. The YAG, Argon & CO2 lasers are expensive compared to Diode technology. QuickLase have developed the use of efficient, reliable diodes in our systems, and as we are the manufacturer by buying direct you "lase out" the middleman for the best value on the market.

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
Troughing / Toma de impresiones

QuickLase® Dual 810+980 nm

The QuickLase Dual system is the only system in the world that combines ultimate ablation and excellent coagulation. The QuickLase Dual combines the wavelength of 810 nanometre, (provides better coagulation and would be preferred in larger area surgeries involving more vascularised tissue), and 980 nanometre (provides better ablation and narrower area involvement and would be preferred in areas where the tissue is less vascularised and/or in areas where narrower tissue involvement is indicated).

The histological evaluation 810nm vs. 980nm wavelength laser radiation presented to the US FDA concluded that mixing the two wavelengths provides definite advantages. General recommendations can be advanced to practitioners on the preferred mixture of the two wavelengths, based on tissue types and dependent procedures.

In figure 2, the vertical axis represents diameter in millimetres. The horizontal axis represents the mixture of wavelengths starting on the left with 100% 810nm and 0% 980nm, then to the right of the graphic with 0% 810nm and 100% 980nm. The top line is the coagulation zone and the bottom line is the ablation zone. The most astounding finding is that when the wavelengths are mixed, refer to the middle portion of both lines, that the effects are greatly magnified: a mixture of the two wavelengths provide superior coagulation and ablation effects when compared to a single wavelength alone.

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The QuickLase Dual is the only system in the market with a FiberCaddy™.

The FiberCaddy™ system is sterilisable complete with both hand piece and fibre optic cable at the same time.

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The beauty of dual wavelengths

Diode lasers are compact, competitively priced, versatile additions to the dentists repertoire predominantly used in performing soft tissue management applications.

The biggest advantage of the diode laser is that when used at the correct wavelength it has a high level of absorption by the haemoglobin, resulting in superb haemostasis, an easy to see bloodless work field and reduced post operative pain.


Figure 1 contains a triangular image that represents the scope and intensity of the coagulation effect of the various wavelength (B), a triangular image that represents the scope and intensity of the ablation effect of the various wavelength (C), and the scale which defines the wavelength measurement in nanometre and the laser medium responsible for generating that particular wavelength (D). As one can see:

- 1.- The shorter wavelengths have a greater effect on coagulation. This is because the shorter wavelengths affect coagulation factors in the blood, directly.
- 2.- The longer wavelengths have a greater effect on tissue removal or ablation. This is because the longer wavelengths are absorbed by water.
- 3.- The on-tissue effects of these wavelengths have been known for decades. Carbon Dioxide (CO₂) lasers produce some of the fastest and most complete ablation known, and Argon lasers produce some of the fastest and most complete photocoagulation witnessed. However, Carbon Dioxide doesn't coagulate and neither do Argon lasers ablate. Before the QuickLase Dual, only the Nd:YAG lasers with frequency doubling, supplying both 532nm and 1062nm wavelengths could provide a complete package of photocoagulation and ablation. The problem with this type of system is the price, costing tens of thousands.
- 4.- As one re-examines figure 1, it is noted that two relatively inexpensive and easy to package alternatives exist; these are the 810nm and 980nm Diode lasers. Many current laser systems output either 810 or 980nm, 810 for better coagulation or 980 for better ablation.

The Quicklase Dual is the only Soft-Tissue Laser System on the market that combines both wavelengths, providing incredible results at an unbelievably low price.

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Advantages

- Minimal or no anaesthetic
- Bloodless
- Minimal suturing
- No collateral damage to the soft tissue
- No damage to the dental pulp
- Can be utilised around implants
- Antibacterial
- Eliminates endotoxin
- Positive treatment for periodontal disease
- Highly regarded technology by the patient.
- No damage to hard tissue

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